



Flask Scrubber® & Steam Scrubber®

实验室玻璃器具清洗机

通用指标

NEW•水温达到 80°C (115V) 93°C(230V)清洗并快速干燥, 空气干燥系统, 快速烘干。

NEW•LCD显示屏可以看到用户设置的循环过程, 并监视整个操作过程(包括内部水温)。

NEW•水循环速率: 363升/分钟 (115V/60Hz和230V/50Hz) 的循环用水。

NEW•双泵设计: 一个清洗, 一个排放, 以减少交叉污染的可能性。

■ 内置式纯净水泵把非加压或加压的纯净水注入6个纯净水喷嘴。

■ 蒸汽发生器对残留物进行熏蒸、软化以清洗玻璃器皿。

- 洗涤剂盒, 用来存放粉末状洗涤剂, 或液体洗涤剂。
- 清洗盒用来存放清洗过程中所需的清洗剂。
- 内置强力空气干燥时间到99分。

NEW•工厂设置的循环程序, 115V 的含 7 个, 230V 增加 3 个, 一共 10 个

- 2 个用户设定的循环程序。

NEW•用户设定延时启动 (到8小时)。

- 304 号抛光不锈钢门, 箱体, 上下清洗臂。
- 一整块不锈钢网孔过滤器, 过滤残留物。

NEW•中央无塔架, 容量增大

- 壹年质保期
- 美国制造
- 在线提供 IQ/OQ 文件, 联系厂家确认服务

长颈瓶清洗机独有特点

NEW■ 304 号不锈钢 36 头支架, (包括 26 个可拆卸的长茎杆清洗长颈瓶和 10 个可拆卸的小支杆清洗小试管及小玻璃容器)

NEW■ 通过支杆进行强力真空干燥

- 提供上支架选件, 含 30 个可拆卸的长茎杆和附加清洗臂
- 提供上下支架可清洗表面皿, 培养管, 其他特殊的玻璃器皿

广口瓶清洗机独有特点

NEW•304 号不锈钢上下支架, 上支架有框

NEW•提供多功能附件---上下多头支架 (没有真空干燥功能)

长颈瓶增强型清洗机独有特点

NEW■ 304 号不锈钢 36 头支架, (包括 26 个可拆卸的长茎杆清洗长颈瓶和 10 个可拆卸的小支杆清洗小试管及小玻璃容器)

- HEPA 不锈钢网孔过滤器经过支杆强力空气干燥
- 水的电导率监控器以确认漂洗进水水质和清洗效果
- RS-232 接口可连接用户端的计算机或记录仪
- 蠕动泵自动液体清洗剂
- 自动酸中和剂

NEW!
Flask Scrubber &
Steam Scrubber racks
now interchangeable



Flask Scrubber Laboratory Glassware Washer 4420430
shown with optional Upper Spindle Rack 4595600



NEW•通过显示器提示信息 and 1 分钟警铃报警

- 透明安全玻璃视窗,
- 内置 25 瓦照明灯, 便于随时观察洗涤过程,
- 可选配带 30 个长茎杆的上支架
- 可选配用于广口瓶清洗的上下支架可清洗培养管, 和其他独特的玻璃器皿

选件和附件

NEW • 柜式和立式 (长颈瓶增强型只有立式的, 但可转换成柜式)

NEW • 移动式转换部件 (长颈瓶机型和广口瓶机型)

NEW • 液体清洁剂分液部件 (长颈瓶增强型清洗机标配)

- 观察视窗, 及内部光源 (长颈瓶增强型清洗机标配)

NEW • 排水系统冷却排水包可清洁外排污水排放达排污法规要求

- 底座

NEW • 多种特殊用途的插件包括16位移液管插件

- Labconco独有特点



Protecting your
laboratory environment
LABCONCO.

 照生有限公司
LIGHT ACE HK LTD.

北京海淀区西三环北路72号, 世纪经贸大厦A座706室, 100037
电话: 010-51798501/02/03/04/05/06 传真: 010-51798500

FlaskScrubber® & SteamScrubber® 实验室玻璃器具清洗机

FlaskScrubber® Laboratory Glassware Washers

长颈瓶清洗机

Catalog #	Style	Electrical Requirements	Internal Water Temperature	With Viewing Window and Light	Overall Dimensions with Door Closed
4420320	Undercounter	115 volts, 60 Hz, 16 Amps	82° C (180° F)	No	24.1" w x 27.4" d x 34.1-36.1" h
4420321	Undercounter	230 volts, 50/60 Hz, 12 Amps	93° C (199° F)	No	24.1" w x 27.4" d x 34.1-36.1" h
4420330	Undercounter	115 volts, 60 Hz, 16 Amps	82° C (180° F)	Yes	24.1" w x 27.4" d x 34.1-36.1" h
4420331	Undercounter	230 volts, 50/60 Hz, 12 Amps	93° C (199° F)	Yes	24.1" w x 27.4" d x 34.1-36.1" h
4420420	Freestanding	115 volts, 60 Hz, 16 Amps	82° C (180° F)	No	24.2" w x 27.5" d x 34.2-36.2" h
4420421	Freestanding	230 volts, 50/60 Hz, 12 Amps	93° C (199° F)	No	24.2" w x 27.5" d x 34.2-36.2" h
4420430	Freestanding	115 volts, 60 Hz, 16 Amps	82° C (180° F)	Yes	24.2" w x 27.5" d x 34.2-36.2" h
4420431	Freestanding	230 volts, 50/60 Hz, 12 Amps	93° C (199° F)	Yes	24.2" w x 27.5" d x 34.2-36.2" h

SteamScrubber® Laboratory Glassware Washers

广口瓶清洗机

Catalog #	Style	Electrical Requirements	Internal Water Temperature	With Viewing Window and Light	Overall Dimensions with Door Closed
4400320	Undercounter	115 volts, 60 Hz, 16 Amps	82° C (180° F)	No	24.1" w x 27.4" d x 34.1-36.1" h
4400321	Undercounter	230 volts, 50/60 Hz, 12 Amps	93° C (199° F)	No	24.1" w x 27.4" d x 34.1-36.1" h
4400330	Undercounter	115 volts, 60 Hz, 16 Amps	82° C (180° F)	Yes	24.1" w x 27.4" d x 34.1-36.1" h
4400331	Undercounter	230 volts, 50/60 Hz, 12 Amps	93° C (199° F)	Yes	24.1" w x 27.4" d x 34.1-36.1" h
4400420	Freestanding	115 volts, 60 Hz, 16 Amps	82° C (180° F)	No	24.2" w x 27.5" d x 34.2-36.2" h
4400421	Freestanding	230 volts, 50/60 Hz, 12 Amps	93° C (199° F)	No	24.2" w x 27.5" d x 34.2-36.2" h
4400430	Freestanding	115 volts, 60 Hz, 16 Amps	82° C (180° F)	Yes	24.2" w x 27.5" d x 34.2-36.2" h
4400431	Freestanding	230 volts, 50/60 Hz, 12 Amps	93° C (199° F)	Yes	24.2" w x 27.5" d x 34.2-36.2" h

FlaskScrubber® Vantage® Series Laboratory Glassware Washer

长颈瓶增强型清洗机

Catalog #	Style	Electrical Requirements	Internal Water Temperature	With Viewing Window and Light	Overall Dimensions with Door Closed
4540031	Freestanding*	230 volts, 50/60 Hz, 12 Amps	93° C (199° F)	Yes	32.2" w x 27.5" d x 34.2-36.2" h

* Convertible to undercounter



SteamScrubber Laboratory Glassware Washer 4400330 includes a viewing window and interior light.



Labconco Corporation
8811 Prospect Avenue
Kansas City, MO 64132-2696
800-821-5525

FAX: 816-363-0130
E-MAIL: labconco@labconco.com
HOME PAGE: www.labconco.com 816-333-8811 or

© 2007 by Labconco Corporation. Printed in the U.S.A. Product design subject to change without notice.



北京海淀区西三环北路72号, 世纪经贸大厦A座706室, 100037
电话: 010-51798501/02/03/04/05/06 传真: 010-51798500